

VistaLab Catalog #	Model	Color	Accuracy*	Precision*
1141C	5µL pipette	Silver	±5.00%	3.74%
1051C	10µL pipette	Orange	±5.00%	2.60%
1142C	15µL pipette	Silver	±5.00%	1.73%
1052C	20µL pipette	Black	±3.00%	1.35%
1053C	25µL pipette	White	±3.00%	1.12%
1143C	30µL pipette	Silver	±3.00%	0.93%
1144C	40µL pipette	Silver	±3.00%	0.70%
1054C	50µL pipette	Green	±2.00%	0.80%
1146C	60µL pipette	Silver	±2.00%	0.67%
1147C	70µL pipette	Silver	±2.00%	0.57%
1148C	75µL pipette	Silver	±2.00%	0.53%
1149C	80µL pipette	Silver	±2.00%	0.50%
1150C	90µL pipette	Silver	±2.00%	0.44%
1055C	100µL pipette	Blue	±1.00%	0.54%
1151C	120µL pipette	Silver	±1.00%	0.45%
1152C	150µL pipette	Silver	±1.00%	0.36%
1056C	200µL pipette	Red	±1.00%	0.30%
1056CLT	200µL pipette <i>(large tip)</i>	Red	±1.00%	0.30%
1057C	250µL pipette	Purple	±1.00%	0.72%
1058C	300µL pipette	Brown	±1.00%	0.60%
1153C	400µL pipette	Silver	±1.00%	0.50%
1154C	450µL pipette	Silver	±1.00%	0.44%
1059C	500µL pipette	Yellow	±1.00%	0.40%
1159C	600µL pipette	Silver	±1.00%	1.00%
1162C	700µL pipette	Silver	±1.00%	0.86%
1157C	750µL pipette	Silver	±1.00%	0.80%
1160C	800µL pipette	Silver	±1.00%	0.75%
1163C	900µL pipette	Silver	±1.00%	0.67%
1060C	1000µL pipette	Gray	±1.00%	0.60%

