

MLA Precision Models

VistaLab Catalog #	Model	Color	Accuracy*	Precision*
1101	5µL pipette	Silver	±5.00%	3.74%
1021	10µL pipette	Orange	±5.00%	2.60%
1102	15µL pipette	Silver	±5.00%	1.73%
1022	20µL pipette	Black	±3.00%	1.35%
1023	25µL pipette	White	±3.00%	1.12%
1103	30µL pipette	Silver	±3.00%	0.93%
1104	40µL pipette	Silver	±3.00%	0.70%
1024	50µL pipette	Green	±2.00%	0.80%
1106	60µL pipette	Silver	±2.00%	0.67%
1107	70µL pipette	Silver	±2.00%	0.57%
1108	75µL pipette	Silver	±2.00%	0.53%
1109	80µL pipette	Silver	±2.00%	0.67%
1110	90µL pipette	Silver	±2.00%	0.57%
1025	100µL pipette	Blue	±1.00%	0.54%
1111	120µL pipette	Silver	±1.00%	0.36%
1112	150µL pipette	Silver	±1.00%	0.36%
1026	200µL pipette	Red	±1.00%	0.30%
1026LT	200µL pipette <i>(large tip)</i>	Red	±1.00%	0.30%
1027	250µL pipette	Purple	±1.00%	0.72%
1028	300µL pipette	Brown	±1.00%	0.60%
1113	400µL pipette	Silver	±1.00%	0.50%
1119	450µL pipette	Silver	±1.00%	0.44%
1029	500µL pipette	Yellow	±1.00%	0.40%
1115	600µL pipette	Silver	±1.00%	1.00%
1505	700µL pipette	Silver	±1.00%	0.86%
1117	800µL pipette	Silver	±1.00%	0.75%
1118	900µL pipette	Silver	±1.00%	0.67%
1030	1000µL pipette	Gray	±1.00%	0.60%

